

XRAYVISION DCV

XrayVision DCV™ and its related tools; components; and plug-ins are designed for 32- or 64-Bit Windows® XP (SP2 or above), Windows Vista®, Windows® 7, Windows Server® 2003, and Windows Server® 2008 operating systems including Terminal Services and Citrix®.

Due to the enterprise nature of XrayVision DCV, server requirements are directly correlated with the size and expected growth of the database and the expected number of users that will be using the application simultaneously. Please contact Apteryx Technical Support for specific system requirements for your site at 877-APTERYX (877-278-3799)

Ideal System Recommendations

- 1.0 GHz or higher processor (faster is always better, especially for servers)
- 256 MB of RAM (again, more is always better)
- 40GB Hard Drive
- Name-brand video capture card* (ATI, Hauppauge)
- CD-R drive
- 100Mbps or higher Local Area Network (TCP/IP preferred)

Minimum System Requirements

Although we don't suggest these specifications for today's digital dental office, Apteryx has set the following as minimums to operate XrayVision DCV:

- 400Mhz or higher processor
- 64MB of RAM (128MB or higher is optimum)
- CD-ROM drive
- Local Area Network for networked installations (TCP/IP preferred; a 100Mbps network is optimum)
- A name-brand video capture card is required for image capture from analog video and intra-oral cameras. Digital sources typically use FireWire or USB for import.



Copyright © 2004-2011 Apteryx, Inc. Apteryx® is a registered trademark of Apteryx, Inc.
All rights reserved. All other trademarks and tradenames are the property of their respective owners.

Controlled Document: Printed or downloaded copies are considered uncontrolled and are for reference only.

REVISION HISTORY

Revision Level	Revision Date	Description of Revision	Effective Date	Revision Author	Approver(s) (Initials)	Doc Change Task#
DRAFT	9/11/2010	Draft	NA	KP	NA	
1.0	01/04/2010	Initial Release		KP		
2.0	10/06/2011	Draft – Doc formatting & Doc name changed to comply with new naming convention.	10/11/2011	SN	PW	85

*** END OF DOCUMENT ***