

## **XRAYVISION DCV**

---

XrayVision DCV and its related tools components, and plug-ins are designed for 32- and 64-Bit Windows® 8, Windows® 10, Windows Server® 2008 / 2008 R2, Windows Server® 2012, and Windows Server® 2016 operating systems including Remote Desktop Services (formerly Terminal Services) and Citrix®

Due to the enterprise nature of XrayVision DCV, server requirements are directly correlated with the size and expected growth of the database and the expected number of users that will be using the application simultaneously. Please contact Apteryx Technical Support for specific system requirements for your site at 877-APTERYX (877-278-3799).

### **Ideal System Recommendations**

---

- 2.0 GHz or higher processor (faster is always better, especially for servers)
- 4GB of RAM (again, more is always better)
- 500 GB Hard Drive
- CD-R/DVDR drive
- 100/1000Mbps or higher Local Area Network (TCP/IP preferred)

### **Minimum System Requirements**

---

The system must conform to the minimum requirements specified by Microsoft® for the operating system installed.



Copyright © 2004-2020 Apteryx, Inc. Apteryx® is a registered trademark of Apteryx, Inc.  
All rights reserved. All other trademarks and tradenames are the property of their respective owners.

**Controlled Document:** Printed or downloaded copies are considered uncontrolled and are for reference only.

## REVISION HISTORY

Revision Level	Revision Date	Description of Revision	Effective Date	Revision Author Initials	Approver Initials	Doc Change Task#
1.0		Initial Release				
2.0		Draft – Doc formatting & Doc name changed to comply with new naming convention.	10/11/2011			85
2.1		Under "Ideal System Requirements" the following has changed: <ul style="list-style-type: none"> <li>- 1.0 GHz is changed to 2.0GHz</li> <li>- 256MB of RAM is changed to 4GB of RAM</li> <li>- 40GB Hard Drive is changed to 500GB</li> <li>- Name-brand video capture card has been deleted</li> <li>- CD-R is changed to CD-R/DVDR</li> <li>- 100Mbps is changed to 100/1000Mbps</li> </ul> "Minimum System Requirements" has been completely rewritten.	08/04/2014			166
2.2	04/05/2018	Removed Windows XP from "designed for" statement, added Windows Server 2016.	04/16/2018	MG	PW	788
2.3	01/27/2020	Removed Windows 7 from "designed for" statement.	01/27/2020	MG	PW	884
2.4	01/27/2020	Removed Windows Vista from "designed for" statement.	01/27/2020	MG	PW	884